



**PRODUCT NO.: 8760**



## **POWDER DETERGENT (Labconco)**

**Non-Foaming, 140 Loads, Contains Phosphates, For Labconco SteamScrubbers**

**DESCRIPTION:** LabSolutions Non-Foaming Powder Detergent is specially formulated for Labconco glassware washers. Removes lab contaminants such as grease, agar, blood and protein digestates. Rinses completely. Non-foaming formula does not promote leaks. A 1% solution has a pH of approximately 12.5. Contains phosphates. Non-hazardous formula may be shipped anywhere. Approximately 140 loads. LabSolutions Powder Detergent, Low-Foaming Liquid Detergent, and Neutralizing Acid Rinse are specially formulated for Labconco SteamScrubbers and FlaskScrubbers. Each removes a wide range of laboratory contaminants and rinses completely.

### **TECHNICAL SPECIFICATIONS:**

- Weight: 10.0 lbs
- Weight metric: 4.5 kg
- Number of Loads: 140

### **ADDITIONAL INFORMATION:**

#### **Applicable Products:**

- 7496-01: Freestanding SteamScrubber
- 8654-01: Undercounter SteamScrubber