



**PRODUCT NO.: 8762**



### **MULTI-PIN INSERT (Labconco)**

7.8"W x 19.3"D x 5.2"H, 48 Pins, Stainless Steel, For Labconco SteamScrubbers

**DESCRIPTION:** Forty-eight pins accommodate a variety of labware. Type 304 stainless steel. Stainless steel pins have plastic tips to protect glassware from scratching. Stainless steel Basket Inserts fit in Standard Open Racks and Spindle Racks.

#### **TECHNICAL SPECIFICATIONS:**

- Weight: 2.0 lbs
- Weight metric: 0.9 kg
- Dimensions, W x D x H (in / cm): 7.8 x 19.3 x 5.2 / 20 x 49 x 13
- Applicable Rack Types: Spindle, Standard
- Rack Location: Lower, Upper

#### **ADDITIONAL INFORMATION:**

##### **Applicable Products:**

- 7496-01: Freestanding SteamScrubber
- 8654-01: Undercounter SteamScrubber