



**PRODUCT NO.: 8980**



## **EZ-DECAPPER, WHEATON (13mm, Green Grip)**

For Standard Seals

**DESCRIPTION:** Standard seals include Center Tear-Out, Open Top, Center Tear-Off. For Flip-Off Seals, once the Flip Top Seal is removed, the seal becomes a Standard Open-Top Seal.

Wheaton E-Z Decappers feature cushioned ergonomic handles that aid in reducing hand fatigue and adds a high degree of comfort for the user. Each decapper size is labeled for quick identification. Other features include polished decapping jaws for consistent seal removal, a short handle and redesigned spring assembly that further aids in reducing hand fatigue. E-Z Decappers can be autoclaved at 121°C. Autoclaving may change appearance and or reduce product life.

### **FEATURES:**

- Used to remove aluminum seals to bottles and vials with a crimp / serum finish
- Cushioned ergonomic handle reduces hand fatigue
- Labeled for quick size identification
- Polished decapping jaws provide consistent seal removal
- Autoclavable at 121°C

### **ADDITIONAL INFORMATION:**

**Warranty:** 2 years