



PRODUCT NO.: 8981



EZ-CRIMPER, WHEATON (13mm, Black)

For Flip-Off Seals, No Grip

DESCRIPTION: These crimpers are used for Flip-Off or Flip Cap Seals only.

FEATURES:

- Used to attach aluminum Flip-Off seals to bottles and vials with a crimp/serum finish
- Polished crimping jaws for consistent sealing
- Spring assembly that further aids in reducing hand fatigue

ADDITIONAL INFORMATION:

Warranty: 1 year