



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 8994



Vacuum Pump, 195L/Min 115V (Freezone Sys)

DESCRIPTION: This vacuum pump provides the deep vacuum needed for good freeze-drying results.

FEATURES:

- ✔ Designed for use with aqueous samples
- ✔ Two position mode selector: High vacuum and high throughput
- ✔ Gas ballast with three positions: closed, low flow and high flow
- ✔ Totally enclosed and fan cooled single phase direct drive motor
- ✔ Isolation valve seals the inlet to prevent oil and air contamination of the system in the event of power failure
- ✔ Comes with four 1 L Vacuum pump oil bottles and one exhaust filter with oil mist and odor filter elements

TECHNICAL SPECIFICATIONS:

Weight (lbs/kg)	65/29.5
Dimensions, W x D x H (in/cm)	6.5 x 18.4 x 10.4/15.8 x 47.0 x 26.1
Electrical	115V, 50/60Hz, 4.6 amps
Sample Contents	Water only
Pump Type	Rotary Vane Vacuum Pump
Works w/ These Products	FreeZone 6 Liter Freeze Dryers
Displacement Capacity at 60Hz (50 Hz)	195 (163) liters/mon
Ultimate Vacuum	2×10^{-3} mBar

APPLICABLE PRODUCTS:

- ✔ 8995: FreeZone Stoppering Tray Dryer 115V (**REQUIRED**)
- ✔ 8996: FreeZone Plus Freeze Dry Sys. 6L Cascade (**REQUIRED**)

WARRANTY: 1 year

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM