



**PRODUCT NO.: 9061**








### **EZ- CRIMPER (Wheaton)**

13mm, For Standard Seals\*, Blue Grip

**DESCRIPTION:** Wheaton E-Z Crimpers feature cushioned ergonomic handles that aid in reducing hand fatigue and adds a high degree of comfort for the user. Each crimper size is labeled for quick identification. Other features include polished crimping jaws for consistent sealing, a short handle and redesigned spring assembly that further aids in reducing hand fatigue. E-Z Crimpers can be autoclaved at 121°C. Autoclaving may change appearance and or reduce product life.

\*Standard seals include Center Tear-Out, Open Top, Center Tear-Off. This does not include Flip Off or Flip Cap seals. This specific crimper will not work for those seals.

#### **FEATURES:**

-  Used to remove aluminum seals to bottles and vials with a crimp / serum finish
-  Cushioned ergonomic handle reduces hand fatigue
-  Labeled for quick size identification
-  Polished crimping jaws provide consistent sealing
-  Autoclavable at 121°C

**Warranty:** 2 years