



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 9170




Vials, Clear Glass, 5mL (Sterile) 20mm

DESCRIPTION: All vials are sterile and pyrogen-free. Prior to vial assembly, vials are dry-heat sterilized and stoppers are steam-sterilized. The vials/stoppers are assembled on an automated line located in a Class 100 suite. The vials are subsequently capped in a Class 100,000 environment. Lot release sterility testing conforms to requirements set forth in USP.

TECHNICAL SPECIFICATIONS:

Filling Volume (mL)	5.0
Height (mm)	36 ± 0.5
Thickness of bottle-lip (mm)	3.6 ± 0.3
Height of bottle-neck (mm)	4.0 ± 0.3
Effective height of bottle-body (mm)	25.0 ± 1.0
Outer diameter of bottle-mouth (mm)	19.8 ± 0.2
Inner diameter of bottle-mouth (mm)	12.6 ± 0.2
Outer diameter of bottle-neck (mm)	15.8 ± 0.3
Body diameter (mm)	22.0 ± 0.2
Material	Borosilicate glass
Cap Colour	Silver
Stopper	Synthetic Chlorobutly, latex-free
Seal	Open-top



USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM