



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 9200



CIMAREC+ STIRRING HOT PLATE, THERMO SCIENTIFIC

10.5" x 10.5", Digital

DESCRIPTION: This digital stirring hotplate offers precise stirring control, safety and temperature performance for routine protocols.

FEATURES:

- ✔ Slow-speed stirring, consistent speed control, and strong magnetic coupling
- ✔ HOT TOP warning system to protect from accidental burns with prominent display when heating surface is above 50 °C (122 °F)
- ✔ Automatic over-temperature protection
- ✔ Easy-to-read LED display for heating and stirring
- ✔ Raised display design protects electronics from spills
- ✔ Stir bar included

WARRANTY: 3 years

TECHNICAL SPECIFICATIONS:

External Dimensions, W x D x H (cm / in)	28.8 x 43.8 x 10.5 / 11.3 x 17.2 x 4.1
Top Plate Dimensions, W x D (cm / in)	26.7 x 26.7 / 10.5 x 10.5
Weight (Kg / lb)	5.9 / 13.0
Top Plate Material	Ceramic
Temperature Range (°C / °F)	30 – 400 (86 – 752)
Stirring Range (RPM)	50 – 1500
Conformity	EN 61000-3-2 EN 61000-3-3 EN 61326-1 EN 61010-1 EN 61010-2-010 EN 61010-2-051