



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 9240**




## Serum Bottle, Clear Glass 50mL, 20mm (USP Type 1)

**DESCRIPTION:** This serum bottle is ideal for long and short-term sample storage, lyophilization and vaccine / injectable drug containers. By matching the tube neck size to the corresponding seal and stopper, a crimper can be used to seal the serum accordingly. These clear bottles are made of borosilicate glass conform to USP Type I requirements. Shrink-wrapped partitioned modules reduce breakage.

### TECHNICAL SPECIFICATIONS:

Filling Volume (mL)	50
Height, H (mm)	70 +0.70/-0.80
Outer diameter of body, D1 (mm)	42.5 +0.70/-0.80
Outer diameter of lip, D2 (mm)	20 0.10/-0.30
Inner diameter of lip, D3 (mm)	12.6 ± 0.20
Outer diameter of neck, D4 (mm)	17
Inner diameter of neck, D5 (mm)	11.5
Material	Borosilicate glass



**USA | CANADA | AUSTRALIA PACIFIC | EUROPE**

[WWW.MEDISCA.COM](http://WWW.MEDISCA.COM)