



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 9497




Clear Glass Vials (Sterile), 2mL, 13mm

DESCRIPTION: This serum bottle is ideal for long and short-term sample storage, lyophilization and vaccine / injectable drug containers. By matching the tube neck size to the corresponding seal and stopper, a crimper can be used to seal the serum accordingly. These clear bottles are made of borosilicate glass conform to USP Type I requirements. Shrink-wrapped partitioned modules reduce breakage.

TECHNICAL SPECIFICATIONS:

Filling Volume (mL)	2
Height (mm)	35 ± 0.5
Thickness of bottle-lip (mm)	3.6 ± 0.2
Height of bottle-neck (mm)	3.5 ± 0.3
Effective height of bottle-body (mm)	23 ± 1.0
Outer diameter of bottle-mouth (mm)	12.8 ± 0.2
Inner diameter of bottle-mouth (mm)	6.9 ± 0.2
Outer diameter of bottle-neck (mm)	10.0 ± 0.3
Body Diameter (mm)	16 ± 0.15
Material	Borosilicate glass
Cap Colour	Silver
Stopper	Synthetic Chlorobutly, latex-free
Seal	Open-top



USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM