



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING



PRODUCT NO.: 9577



CONICAL CYLINDER, GRADUATED PLASTIC (250 mL)
Graduation Range: 50-250 mL, Dual Scale, Clear w/Black Base

DESCRIPTION: This conical cylinder provides excellent chemical resistance to its content. It is graduated in ounces and milliliters with no meniscus to confuse readings. It conforms to the volume and accuracy requirements of ASTM E1094.

FEATURES:

-  Autoclavable (will affect accuracy)
-  Graduated

TECHNICAL SPECIFICATIONS:

Capacity (mL / oz)	250 / 8
Type	Plastic Graduated Cylinder
Limit of Error (mL)	±3.5
Material	Polymethylpentane (PMP)
Compliance	ASTM E1094