



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 9619






WEIGH PAPER

W 4" × D 4"

DESCRIPTION: This weigh paper protects the pans of weighing balances from damage and provides a convenient way to transfer weighed materials to another container.

New weighing paper should be folded along its diagonals, flattened and then placed on the pan for different drugs to prevent contamination.

FEATURES

-  Nitrogen-free
-  Smooth surface, allowing solids to flow freely
-  Supports relatively heavy weight