



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 9709






Syringes, Luer Lock, 50mL (Sterile)

DESCRIPTION: This luer-lock tip syringe features a clear barrel with bold scale marking, tapered plunger rod for ease of aspiration, positive plunger rod-stop and an added BD Luer-Lock for increased secure connection.

It does not come with a needle.

FEATURES:

-  Latex-free
-  Disposable
-  Ethylene Oxide (EtO) sterilized

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM