



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 9710






LUER-LOK TIP SYRINGE (30 mL, Sterile)

Graduation Indication: 1.0 mL, Clear

DESCRIPTION: This luer-lock tip syringe features a clear barrel with bold scale marking, tapered plunger rod for ease of aspiration, positive plunger rod-stop, and an added BD Luer-Lok thread for increased secure connection.

It does not come with a needle.

FEATURES:

-  Latex-Free
-  Disposable
-  Ethylene Oxide (EtO) sterilized