



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 9786**



## **Mortar & Pestle, 473.2mL (16oz)**

Heavy Clear Flint Glass

**DESCRIPTION:** This mortar and pestle have been used for many years, in the pharmaceutical compounding of drugs, to grind crystals and large particles into fine powders. The mortar and pestle are both made of heavy clear flint glass which provides a clear view of its contents. The mortar has a heavy molded base and lip for steady support and easy pouring.

### **TECHNICAL SPECIFICATIONS:**

Capacity (mL/oz)	473.2 / 16
Dimensions, H x ID (mm)	135 x 120
Material	Heavy Glass*

\* Please note that this product is not considered dishwasher safe.

**USA | CANADA | AUSTRALIA PACIFIC | EUROPE**

[WWW.MEDISCA.COM](http://WWW.MEDISCA.COM)