



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 9790



Mortar & Pestle, 29.5mL (1oz)

Heavy Clear Flint Glass

DESCRIPTION: This mortar and pestle have been used for many years, in the pharmaceutical compounding of drugs, to grind crystals and large particles into fine powders. The mortar and pestle are both made of heavy clear flint glass which provides a clear view of its contents. The mortar has a heavy molded base and lip for steady support and easy pouring.

TECHNICAL SPECIFICATIONS:

Capacity (mL/oz)	29.5 / 1
Dimensions, H x ID (mm)	70 x 60
Material	Heavy Glass*

* Please note that this product is not considered dishwasher safe.

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM