



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 9793



MORTAR & PESTLE (145 mL / 4.9 oz)

Porcelain

DESCRIPTION: This mortar and pestle have been used for many years, in the pharmaceutical compounding of drugs, to grind crystals and large particles into fine powders. The mortar and pestle are both made of porcelain.

The outside and rim of the mortar and the pestle handle are glazed for the purpose of dirt and stain resistance.

TECHNICAL SPECIFICATIONS:

Mortar Diameter (mm)	100
Capacity (mL / oz)	145 / 4.9
Material	Porcelain