




**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING




**PRODUCT NO.: 9802**



## **Supposistrip Mold (12 × 1.76 mL Cavities, Pink)**

**DESCRIPTION:** This re-usable, flexible supposistrip mold is designed to produce suppositories that can be pushed or pinched out of each cavity. The strips can also be trimmed with a pair of scissors, depending on the number of cavities required for the prescription.

### **FEATURES:**

-  Made of Dynaflex, a thermoplastic elastomer designed to have a rubbery feel and soft touch
-  Latex free
-  BPA free

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

[WWW.MEDISCA.COM](http://WWW.MEDISCA.COM)