



PRODUCT NO.: 9831



ERLENMEYER FLASK, NARROW MOUTH (1000 mL)

Graduation Range: 250 – 1000 mL, Graduation Interval: 50 mL

DESCRIPTION: This flask is designed with heavy duty rims to reduce chipping. The uniform wall thickness provides the proper balance between mechanical strength and thermal shock resistance. Approximate graduations are in durable white enamel and an extra-large marking space is provided.

It can be used for several laboratory applications for holding and storage purposes.

TECHNICAL SPECIFICATIONS:

Capacity (mL)	1000
Flash Style	Erlenmeyer
Neck Style	Narrow Mouth
Top Style	Heavy Duty Rim
Rubber Stopper No.	9
Graduation Range (mL)	250 – 1000
Graduation Interval (mL)	50
Height (mm)	~216
Outer Diameter (mm)	~129