



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 9834**



## **ERLENMEYER FLASK, NARROW MOUTH (125 mL)**

Graduation Range: 50 – 125 mL, Graduation Interval: 25 mL

**DESCRIPTION:** This flask is designed with heavy duty rims to reduce chipping. The uniform wall thickness provides the proper balance between mechanical strength and thermal shock resistance. Approximate graduations are in durable white enamel and an extra-large marking space is provided.

It can be used for several laboratory applications for holding and storage purposes.

### **TECHNICAL SPECIFICATIONS:**

Capacity (mL)	125
Flash Style	Erlenmeyer
Neck Style	Narrow Mouth
Top Style	Heavy Duty Rim
Rubber Stopper No.	5
Graduation Range (mL)	50 – 125
Graduation Interval (mL)	25
Height (mm)	~114
Outer Diameter (mm)	~67