



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 9860







CONICAL CYLINDER, GRADUATED GLASS TYPE I (500 mL)

Graduated Range: 100 – 500 mL, Graduation Interval: 25 mL, Pharmaceutical Grade, Dual Scale

DESCRIPTION: This conical cylinder is used to perform all sort of measurements.

FEATURES:

-  Scales in both metric and U.S. customary units
-  "TD" appears on the cylinder and indicates it is calibrated to deliver
-  Designed from ASTM Specification E1094, Type 2 requirements (also NBS Handbook 44)
-  Manufactured from 33-expansion, low extractable borosilicate glass conforming to USP Type I and ASTM E438, Type I, Class A requirements

TECHNICAL SPECIFICATIONS:

Capacity (mL)	500
Height (mm)	243
Graduation Range (mL)	100 – 500
Graduation Interval (mL)	25
Tolerance (mL)	Up to 100 ± 2.2, above 100 – 150 ± 2.8, above 150 – 275 ± 3.4, above 275 – 425 ± 4.1, above 425 ± 4.8