



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 9861



CONICAL CYLINDER, GRADUATED GLASS TYPE I (250 mL)

Graduated Range: 50 – 250 mL, Graduation Interval: 25 mL, Pharmaceutical Grade, Dual Scale

DESCRIPTION: This conical cylinder is used to perform all sort of measurements.

FEATURES:

- ▤ Scales in both metric and U.S. customary units
- ▤ "TD" appears on the cylinder and indicates it is calibrated to deliver
- ▤ Designed from ASTM Specification E1094, Type 2 requirements (also NBS Handbook 44)
- ▤ Manufactured from 33-expansion, low extractable borosilicate glass conforming to USP Type I and ASTM E438, Type I, Class A requirements

TECHNICAL SPECIFICATIONS:

Capacity (mL)	250
Height (mm)	193
Graduation Range (mL)	50 – 250
Graduation Interval (mL)	25
Tolerance (mL)	Up to 50 ± 1.4, above 50 – 100 ± 1.8, above 100 – 175 ± 2.2, above 175 ± 2.8