



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 9864







CONICAL CYLINDER, GRADUATED GLASS TYPE I (25 mL)

Graduated Range: 5 – 25 mL, Graduation Interval: 5 mL, Pharmaceutical Grade, Dual Scale

DESCRIPTION: This conical cylinder is used to perform all sort of measurements.

FEATURES:

-  Scales in both metric and U.S. customary units
-  "TD" appears on the cylinder and indicates it is calibrated to deliver
-  Designed from ASTM Specification E1094, Type 2 requirements (also NBS Handbook 44)
-  Manufactured from 33-expansion, low extractable borosilicate glass conforming to USP Type I and ASTM E438, Type I, Class A requirements

TECHNICAL SPECIFICATIONS:

Capacity (mL)	25
Height (mm)	113
Graduation Range (mL)	5 – 25
Graduation Interval (mL)	5
Tolerance (mL)	@ 5 ± 0.2, above 5 – 15 ± 0.4, above 15 ± 0.6