



PRODUCT NO.: 9865



CONICAL CYLINDER, GRADUATED GLASS TYPE I (10 mL)

Graduated Range: 0.5 – 10 mL, Graduation Interval: 0.5 mL (From 0.5 – 2 mL) & 1 mL (From 2 – 10 mL), Pharmaceutical Grade, Dual Scale

DESCRIPTION: This conical cylinder is used to perform all sort of measurements.

FEATURES:

- ▣ Scales in both metric and U.S. customary units
- ▣ "TD" appears on the cylinder and indicates it is calibrated to deliver
- ▣ Manufactured from 33-expansion, low extractable borosilicate glass conforming to USP Type I and ASTM E438, Type I, Class A requirements

TECHNICAL SPECIFICATIONS:

Capacity (mL)	10
Height (mm)	82
Graduation Range (mL)	0.5 – 10
Graduation Intervals (mL)	0.5 – 2 in 0.5, 2 – 10 in 1
Tolerance (mL)	Up to 3 ± 0.1, above 3 ± 0.2

ADDITIONAL INFORMATION:

- ▣ Does not meet ASTM E1094 nor NBS Handbook 44 requirements and therefore may not be used for prescription work in states which adopt these specifications