



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 9938



BEAKER, GRADUATED GLASS (1000 mL)

Graduation Range: 50–1000 mL, Graduation Interval: 50 mL, Dual Scale, Heavy Duty

DESCRIPTION: This heavy duty beaker is designed to offer good mechanical strength under harsh conditions, such as mechanized washing operations. It has an extra-large marking spot and conforms to ASTM E-960 standards.

TECHNICAL SPECIFICATIONS:

Capacity (mL)	1000
Graduation Range (mL)	50 – 1000
Graduation Interval (mL)	50
Accuracy	± 5% of full capacity
Approximate O.D. × H (mm)	108 × 156