



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING




PRODUCT NO.: 9950



WIPE, KIMTEX (8 3/4" × 16 3/4")

DESCRIPTION: This wipe is ideal for tasks where precision and reliability is essential and can be used to clean and disinfect surfaces in healthcare facilities. It is designed for wiping tasks where low lint or use of chemicals is required. It is resistant to acids, alkalis, solvents, inks and paints. Polypropylene construction is perfect for use with lubricants, oils and greases.

FEATURES:

-  Reusable, heavy-duty wipers for use with oils and solvents
-  Polypropylene construction for cloth-like feel and long-lasting performance
-  Low-lint cleaning performance absorbs quickly for big cleanups

TECHNICAL SPECIFICATIONS:

Dimension	914 x 676 x 889/36 x 27 x 35
Material	Polypropylene
Color	Blue
Delivery system	Pop-Up Box

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM