



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 9973



BEAKER, GRADUATED GLASS TYPE I (2000 mL)

Graduation Range: 200–1800 mL, Graduation Interval: 100 mL, Dual Scale

DESCRIPTION: This beaker meets ASTM E-960 Type I Specifications for classification, design, capacity & dimensions, and markings. This beaker is used to hold and sometimes measure liquids. It is also used to stir, mix and heat liquids and it is made from high-quality borosilicate glass with a low coefficient of expansion.

FEATURES:

- ✔ Easy-pour sprout
- ✔ White, permanent graduations
- ✔ Large marking squares
- ✔ Manufactured with uniform wall thickness resulting in the optimum balance between thermal shock resistance mechanical strength

TECHNICAL SPECIFICATIONS:

Capacity (mL)	1000
Material	High Quality Borosilicate Glass
Graduation Range (mL)	200 – 1800
Graduation Interval (mL)	100
Accuracy	± 5% of full capacity
Approximate O.D. x H (mm)	130 x 188