



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 9977**







### **Beaker, Graduated Glass Type 1 (400mL)**

**Graduation Range: 25-325mL, Graduation Interval: 25mL, Dual Scale**

**DESCRIPTION:** This beaker meets ASTM E-960 Type I Specifications for classification, design, capacity & dimensions, and markings. This beaker is used to hold and sometimes measure liquids. It is also used to stir, mix and heat liquids and it is made from high-quality borosilicate glass with a low coefficient of expansion.

#### **FEATURES:**

-  Easy-pour sprout
-  White, permanent graduations
-  Large marking squares
-  Manufactures with uniform wall thickness resulting in the optimum balance thermal shock resistance mechanical strength

#### **TECHNICAL SPECIFICATIONS:**

Capacity (mL)	400
Material	High Quality Borosilicate Glass
Graduation Range (mL)	25 - 325
Graduation Interval (mL)	25
Accuracy	± 5% of full capacity
Approximate O.D. x H (mm)	78 x 112

**USA | CANADA | AUSTRALIA PACIFIC | EUROPE**

[WWW.MEDISCA.COM](http://WWW.MEDISCA.COM)