



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 9980



BEAKER, GRADUATED GLASS TYPE I (100 mL)

Graduation Range: 20–80 mL, Graduation Interval: 10 mL, Dual Scale

DESCRIPTION: This beaker meets ASTM E-960 Type I Specifications for classification, design, capacity & dimensions, and markings. This beaker is used to hold and sometimes measure liquids. It is also used to stir, mix and heat liquids and it is made from high-quality borosilicate glass with a low coefficient of expansion.

FEATURES:

- ✔ Easy-pour sprout
- ✔ White, permanent graduations
- ✔ Large marking squares
- ✔ Manufactured with uniform wall thickness resulting in the optimum balance between thermal shock resistance mechanical strength

TECHNICAL SPECIFICATIONS:

Capacity (mL)	100
Material	High Quality Borosilicate Glass
Graduation Range (mL)	20 – 80
Graduation Interval (mL)	10
Approximate O.D. × H (mm)	50 × 75