




MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING




PRODUCT NO.: 9999



Supposistrip Mold (12 × 1.35 mL Cavities, Blue)

DESCRIPTION: This re-usable, flexible supposistrip mold is designed to produce suppositories that can be pushed or pinched out of each cavity. The strips can also be trimmed with a pair of scissors, depending on the number of cavities required for the prescription.

FEATURES:

-  Made of Dynaflex, a thermoplastic elastomer designed to have a rubbery feel and soft touch
-  Latex free
-  BPA free

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM