


PRODUCT NO.: 5844



Bottle, Amber Glass, Boston Round, 8oz, 24-400 w/black Phenolic Cap

Description: Amber Boston Round Bottles protect contents from UV rays and are ideal for light sensitive products. The general use bottles are perfect for liquids.

Technical Specifications:

Capacity	8 oz.
Neck finish	24-400
Style	Boston round bottle
Color	Amber
Width (mm)	61.5
Height (mm)	136.9
Glass type	Soda-lime silica

USA | Canada | Australia | International

medisca.com