




MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 6282



CANOPY CONNECTION, LABCONCO (4 ft) For RXPert Double Filtered Balance Systems

DESCRIPTION: This canopy connection allows the enclosure to be thimble-ducted to the outside for additional protection. The canopy's air gap provides a means of balancing air by allowing air to enter the ductwork and mix with the exhausted filtered air from the enclosure.

TECHNICAL SPECIFICATIONS:

Enclosure Width (ft)	4
Duct Dimension, ID × OD (in)	6.1 × 6.2
Shipping Weight (kg / lb)	10 / 21
Material	Powder-coated Steel