



**PRODUCT NO.: 7339**



## **REPLACEMENT BAG FOR BAG-IN/BAG-OUT HEPA/ULPA FILTER** **4ft, For Labconco XPert Filtered Balance Systems**

**DESCRIPTION:** This product is required when changing the HEPA filter on a Labconco XPert Filtered Balance System in order to protect the operator from any hazardous substances that are present in the HEPA filter.

### **FEATURES:**

- Weight: 9.0 lbs
- Weight metric: 4.1 kg
- Dimensions: 96.0" circumference x 96.0" long
- Dimensions metric: 243.8 x 243.8 cm
- Applicable Categories: XPert Balance & Bulk Powder Enclosures, XPert Nano Enclosures
- Applicable Products: XPert Filtered Balance Stations, XPert Filtered Balance Systems, XPert Nano Enclosures
- Enclosure Nominal Width: 4'
- Filter Application: Particulate
- Filter Length (Inches): 90 to 100
- Filter Type: Filter Accessory

### **ADDITIONAL INFORMATION:**

#### **Applicable Products:**

- 8311-01
- 9332-01
- 8797-01
- 8392-01