



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 9857



CYLINDER, GRADUATED GLASS TYPE I, CLASS B (10 mL)

Graduated Range: 0.7 – 10.0 mL, Graduation Interval: 0.1 mL, Pharmaceutical Grade, Double Scale

DESCRIPTION: This graduated cylinder meets ASTM E-1272 Class B, Type I Specifications & USP standards for volumetric glassware. Hexagonal base provides extra stability and all markings are in durable, white enamel. Comes with a protective, polyethylene, hexagonal, bumper guard.

TECHNICAL SPECIFICATIONS:

Capacity (mL)	10
Approximate O.D. x Height (mm)	8 x 182
Graduation Range (mL)	0.7 – 10.0
Graduation Interval (mL)	0.1
Tolerance (mL)	± 0.10