PRODUCT NO.: 8795




Adapters, Press-In Bottle, for 24 mm Bottle

Description: Press-in bottle adapters fit securely in any bottle with the appropriate size opening. The Press-in Bottle Adapters allow a user to easily draw liquid with an oral syringe. The adapters fit flush with the bottle opening so that the original cap can be placed on the bottle.

Features:

- Clear color
- FDA approved plastic material
- Made up of LDPE

Which Press-In adapters will fit my bottle?					
Adapter Stock Code	For Bottle Size	Adapter Specifications			
		Outer Diameter	Inner Diameter	Height	Diagram
8796	20 mm	16.26 mm	11.18 mm	14.78 mm	
8795	24 mm	20.42 mm	15.14 mm	15.41 mm	
8794	28 mm	22.62 mm	17.8 mm	15.03 mm	
6349	33 mm	28.05 mm	22.66 mm	15.38 mm	
Note: MEDISCA cannot guarantee full standardization of all bottles on the market. Please test your adapter for final confirmation of fit.					

Please note that the information contained herein is for informational purposes only. The information is based on recognized compendia, manufacturers' instructions and/or other authoritative sources and texts. MEDISCA takes no responsibility for the validity or accuracy of this information. Customer MUST refer to USP, FDA, Health Canada, TGA and/or their State Board to ensure they comply with the specific regulations and ensure that appropriate procedures are in place.